

Prior to sending a scan, confirm the following:

1. The scan body is correct.
2. The scan body was scanned fully and the bevel is facing buccal or lingual and not interproximal.
3. Contact points are broad rather than points.
4. Insertion direction.

The versatile treatment for your patients.

The NobelReplace implant offers a trusted solution to a wide range of indications, and further flexibility when critical anatomical structures limit space.

Key Indications	Fixation	Material	Workflow	Implant Systems	Implant Connections
Single Unit	Cement or screw retained	Titanium	Digital	3.0/NP/RP/WP	CC/Tri-Ch

To access this program contact information@nutooth.ca for an account username and password or go to www.nutooth.ca and chat online with one of our dental experts.



ANNOUNCEMENT



48 HOURS



Scan received by 2 pm and the case will be shipped within 48 hours.

Now an option on the *nutooth* prescription portal.

nutooth



SEND A CASE

ANTERIOR/POSTERIOR AUTHENTIC PARTS WITH CUSTOM DESIGNS



JUSTIFIED TRUST WITH
THE ORIGINAL IMPLANT-
ABUTMENT CONNECTION

- Perfectly matching components
- Excellent service and support
- Nobel lifetime warranty included

ABUTMENT AND
FCZ MILLED AND
DESIGNED BY

*nu*tooth

NOBEL MILLED
CUSTOM ABUTMENT
AND A FCZ
ANTERIOR/POSTERIOR
SCREW RETAINED
CROWN
ON A DIGITAL SCAN



UNIBASE
ABUTMENTS
AVAILABLE FOR
APPLICABLE
CASES