


ANTERIOR/POSTERIOR AUTHENTIC PARTS WITH CUSTOM DESIGNS



JUSTIFIED TRUST WITH
THE ORIGINAL IMPLANT-
ABUTMENT CONNECTION

- Perfectly matching components
- Excellent Service and support
- Straumann Guarantee*

ABUTMENT AND
FCZ MILLED AND
DESIGNED BY

*nu*tooth

STRAUMANN MILLED
CUSTOM ABUTMENT
AND A FCZ
ANTERIOR/POSTERIOR
SCREW RETAINED
CROWN
ON A DIGITAL SCAN



VARIOBASE
ABUTMENTS
AVAILABLE FOR
APPLICABLE
CASES

Prior to sending a scan, confirm the following:

1. The scan body is correct.
2. The scan body was scanned fully and the bevel is facing buccal or lingual and not interproximal.
3. Contact points are broad rather than points.
4. Insertion direction.

Original components from Straumann make the difference. When you choose original Straumann components, you receive a high-precision product of exceptional quality. The perfect match of Straumann implants and abutments results in excellent fatigue test performance, while the exceedingly thorough inspection process ensures consistent quality across all original Straumann components.

Key Indications	Fixation	Material	Workflow	Implant Systems	Implant Connections
Single Unit	Cement or screw retained	TAN	Digital	BL/BLX/TL	NC/RC/WN/RN RB/WB

To access this program contact information@nutooth.ca for an account user and PW or go to www.nutooth.ca and chat online with one of our dental experts.

*Straumann warranty applies



Straumann® Original
for dentists

When failure is not an option.

ANNOUNCEMENT



48 HOURS




Scan received by 2 pm and the case will be shipped within 48 hours.

Now an option on the *nutooth* prescription portal.

nutooth

 **straumann**



SEND A CASE